COMFORT QUALITY COMPLIANT EFFECTIVE





Fusion[™]

CUSTOM CALIBRATION

- Patient specific calibration
- 3 removable wings for 1mm of protrusion on the bottom
- Fine calibration in .1mm increments with lugs on the top
 Longer range of advancement (8 Emm) leading to forward
- Longer range of advancement (8.5mm) leading to fewer device resets



Classic

ORIGINAL INNOVATION

- Highest quality acrylic
- Medical grade stainless steel ball clasps for retention
- Lingual-less option available for additiona tongue space





VISUAL SIMPLICITY

- E0486 Medicare verified
- Visual calibration indicator makes it easy to measure how far device has been advanced
- SMH BFlex offers superior comfort and retention
- Longer range of advancement (8mm) leading to fewer device resets over the life of the device

SomnoMed® & SomnoDent® are registered trademarks of SomnoMed. Herbst® is a registered trademark of Dentaurum, Inc., Newtown, PA. SM FULL LINE PRODUCT SELL SHEET 903004revC



PREMIUM COMFORT

- Proprietary SMH BFlex inner liner for superior retention
- Drop in fit even for patients with short teeth, little undercut or crown / bridge work



LATEST INNOVATION

- Ease of calibration
- Wide range of calibration (-2mm to +8mm)
- Standard metal-free, also available with ball clasps




POPULAR VALUE

- E0486 Medicare verified
- Standard telescopic mechanism (also available in Shim)
- SMH BFlex offers superior comfort and retention, also available with ball clasps



7460 Warren Parkway, Suite 190, Frisco, TX 75034 🥃 www.somnomed.com

SomnoMed[®] offers a broad range of oral devices for the first line treatment of mild to moderate obstructive sleep apnea.



SomnoDent[®] devices are:

Comfortable for Patients

- 88% of patients reported regular use of their SomnoDent device¹
- Unique SMH BFlex technology provides unparalleled comfort²
- Patients can talk and open/close their mouths and lips while wearing device

High Quality

- FDA 510K cleared and manufactured at an ISO 13485 certified facility
- Customized using the highest quality acrylic which does not discolor or attract odors
- Device life expectancy 3-5 years; and devices have a 3 year warranty against manufacturing defects.³

Uniquely Compliant

- The 1st and only dorsal fin oral device with micro-recorder that objectively records compliance data
- Reports include hours device is worn and patient's supine or non-supine head position

Clinically Effective

- 15+ independent studies demonstrate the device's significant clinical benefits
- 91% of patients reported improvement in sleep quality while wearing SomnoDent¹

¹AM / Respir Crit Care Med Vol 163. Pp 1457-1461, 2001 ²SMH BFlex not available in Classic device ³See website for details

