

Atlantis Compatible Implant Platforms

| Manufacturer | Implant | ATLANTIS Abutment | | ATLANTIS Crown | | ATLANTIS Conus Abutment | |
|---|--|-------------------------|----------|----------------|----------|-------------------------|-------------|
| | | | | | | Custom | Overdenture |
| | | Titanium Gold Shaded Ti | Zirconia | Titanium | Zirconia | Titanium Gold Shaded Ti | Titanium |
| BioHorizons | INTERNAL | | | | | | |
| | Internal 3.5, 4.5, 5.7 | X | X | X | X | X | X |
| | Tapered Internal 3.5, 4.5, 5.7 | X | X | X | X | X | X |
| | EXTERNAL | | | | | | |
| | External Implant 3.5 | X | | | | X | |
| External Implant 4.0, 5.0 | X | X | X | X | X | | |
| Biomet 3i | INTERNAL | | | | | | |
| | Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0 | X | X | X | X | X | X |
| | Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5 | X | X | X | X | X | X |
| | Certain™ XP 4/5, 5/6 | X | X | X | X | X | X |
| | EXTERNAL | | | | | | |
| | MicroMiniplant™ 3.25 | X | | | | X | |
| | Miniplant™ 3.25 | X | X | | X | X | |
| Osseotite™ 3.75, 4.0, 5.0, 6.0 | X | X | | X | X | | |
| Osseotite™ XP 4/5, 5/6 | X | X | | X | X | | |
| Camlog | INTERNAL | | | | | | |
| | Screw-line Implant 3.3 | X | | X | | X | X |
| | Screw-line/Root-line Implant 3.8, 4.3, 5.0, 6.0 | X | X | X | X | X | X |
| DENTSPLY Implants | INTERNAL | | | | | | |
| | ANKYLOS® C/X 3.5, 4.5, 5.5, 7.0 | X | X | | | X | X |
| | OsseoSpeed™ EV 3.0 – Green | X | | X | | X | X |
| | OsseoSpeed™ EV 3.6 – Purple | X | X | X | X | X | X |
| | OsseoSpeed™ EV 4.2 – Yellow | X | X | X | X | X | X |
| | OsseoSpeed™ EV 4.8 – Blue | X | X | X | X | X | X |
| | OsseoSpeed™ EV 5.4 – Brown | X | X | X | X | X | X |
| | OsseoSpeed™/OsseoSpeed™ TX 3.0 – Yellow | X | | X | | X | X |
| | OsseoSpeed™/OsseoSpeed™ TX 3.5, 4.0 – Aqua | X | X | X | X | X | X |
| | OsseoSpeed™/OsseoSpeed™ TX 4.5, 5.0 – Lilac | X | X | X | X | X | X |
| | OsseoSpeed™ TX Profile 4.5, 5.0 | X | | X | | X | X |
| | OsseoSpeed™ 4.0 – Blue (Previous) | X | | | | X | |
| | XiVE® S 3.0 | X | | X | | X | X |
| XiVE® S 3.4, 3.8, 4.5, 5.5 | X | X | X | X | X | X | |
| Keystone Dental | INTERNAL | | | | | | |
| | PrimaConnex™ SD 3.3, 3.5, Genesis™ 3.8 | X | X | X | X | X | X |
| | PrimaConnex™ RD 4.0, 4.1, Genesis™ 4.5 | X | X | X | X | X | X |
| | PrimaConnex™ WD 5.0, Genesis™ 5.5, 6.5 | X | X | X | X | X | X |
| | EXTERNAL | | | | | | |
| | Restore™ SD 3.3, 3.4 | X | | | | X | |
| | Restore™ RD 3.75, 4.0 | X | X | | X | X | |
| Restore™ WD 5.0, 5.5, 6.0 | X | | | | X | | |
| Nobel Biocare | INTERNAL | | | | | | |
| | NobelActive™ NP 3.5 | X | X | X | X | X | X |
| | NobelActive™ RP 4.3, 5.0 | X | X | X | X | X | X |
| | NobelReplace™ NP 3.5 | X | X | X | X | X | X |
| | NobelReplace™ RP 4.3 | X | X | X | X | X | X |
| | NobelReplace™ WP 5.0 | X | X | X | X | X | X |
| | NobelReplace™ 6.0 | X | X | X | X | X | X |
| | NobelReplace™ Conical Connection NP 3.5 | X | X | X | X | X | X |
| | NobelReplace™ Conical Connection NP 4.3, 5.0 | X | X | X | X | X | X |
| | EXTERNAL | | | | | | |
| | Brånemark System™ NP 3.3 | X | | | | X | |
| | Brånemark System™ RP 3.75, 4.0 | X | X | | X | X | |
| | Brånemark System™ WP 5.0, 5.5 | X | X | | X | X | |
| SteriOss Replace™ 3.1, 3.5, 4.3, 5.0, 6.0 | X | | | | X | | |
| SteriOss HL™ 3.25, 3.8, 4.5, 5.0, 6.0 | X | | | | X | | |
| Straumann | INTERNAL | | | | | | |
| | Bone Level 3.3 NC (Narrow CrossFit™) | X | X | X | X | X | X |
| | Bone Level 4.1, 4.8 RC (Regular CrossFit™) | X | X | X | X | X | X |
| | Standard/Standard Plus 4.8 RN (Regular Neck) | X | | X | | X | X |
| | Standard/Standard Plus 4.8 WN (Wide Neck) | X | | X | | X | X |
| | EXTERNAL | | | | | | |
| Standard Plus 3.3 NN (Narrow Neck) | X | | | | X | | |
| Sybron Implant Solutions | EXTERNAL | | | | | | |
| | Endopore 3.5, 4.1, 5.0 - external x 0.050 | X | | | | | |
| | Entegra 3.5, 4.1, 5.0 - external | X | | | | | |
| Zimmer Dental | INTERNAL | | | | | | |
| | Screw Vent™ 3.3, 3.7, 4.5 | X | X | X | X | X | X |
| | Screw-Vent™ 3.3, 3.7, 4.5 friction fit | X | | X | | X | X |
| | Tapered Screw-Vent™ 3.5, 4.5, 5.7 | X | X | X | X | X | X |
| | Tapered Screw-Vent™ 3.5, 4.5, 5.7 friction fit | X | | X | | X | X |
| SwissPlus™ 4.8 | X | | X | | X | X | |

Custom Abutment Compatible Implant Platforms

| | Implant | Zirconia | Titanium |
|---|--|----------|----------|
| BioHorizons Authentic | INTERNAL | | |
| | 3.0, 3.5, 4.5, 5.7 | | X |
| | 3.5, 4.5, 5.7 | X | |
| Nobel Biocare Custom Abutments | INTERNAL | | |
| | Straumann Bone Level - NC, RC | | X |
| | Straumann Synocta - RN, WN | | X |
| | Nobel Active - 3.0, 3.5, 4.3, 5.0 | | X |
| | Nobel Replace - 3.5, 4.3, 5.0, 6.0 | | X |
| | Astra Tech Osseospeed - 3.0, 3.5/4.0, 4.5/5.0 | | X |
| | Biomet 3i Certain - 3.4, 4.0, 5.0, 6.0 | | X |
| | Zimmer Screw Vent - 3.5, 4.5, 5.7 | | X |
| | Straumann Bone Level - NC, RC | X | |
| | Straumann Synocta - RN, WN | X | |
| | Nobel Active - 3.5, 4.3 | X | |
| | Nobel Replace - 3.5, 4.3, 5.0, 6.0 | X | |
| | Astra Tech - 3.0, 3.5/4.0, 4.5/5.0 | X | |
| | EXTERNAL | | |
| | Biomet Certain External Hex - 3.25, 4.0, 5.0, 6.0 | | X |
| Straumann Authentic Custom Abutments | INTERNAL | | |
| | Straumann Bone Level - NC, RC | | X |
| | Straumann Synocta - RN, WN | | X |
| | Straumann Bone Level - NC, RC (one piece abutment) | X | |
| | Straumann Synocta - RN (one piece abutment) | X | |
| | Straumann Bone Level - NC, RC (two piece abutment) | X | |
| Straumann Synocta - RN, WN (two piece abutment) | X | | |

Oral Arts Custom Abutments

| Compatible Implant Platforms | Restorative Platform Sizes | | |
|---|----------------------------|-----------|-------|
| Zimmer | 3.5mm | 4.5mm | 5.7mm |
| BioHorizons | 3.5mm | 4.5mm | 5.7mm |
| Implant Direct Legacy | 3.5mm | 4.5mm | 5.7mm |
| Nobel Active | 3.5mm | 4.3mm | |
| Nobel Replace | 3.5mm | 4.3mm | 5.0mm |
| Biomet 3i Certain | 3.4mm | 4.1mm | 5.0mm |
| Astra | 3.5/4.0mm | 4.5/5.0mm | |
| Straumann Bone Level | NC | RC | |
| <i>3.0 and external hex implants are not available at this time.</i> | | | |